

Anti-Rpl24 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P83731

Item Number: CM0015039



Introduction

Rabbit polyclonal antibody targeting human RPL24 (ribosomal protein L24). Validated for Western blot, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat RPL24. Suitable for translation research.

[Learn More](#)

Parameter	Value
Product Name	RPL24 Rabbit pAb
Remarks/Alias	L24; eL24; HEL-S-310; RPL24
Species	Human
GeneID (Human)	6152
GeneID	6152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-157 of human RPL24 (NP_000977.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P83731
Protein Weight	18kDa
Shipping	Ice bag