


Anti-Rpl26 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P61254

Item Number: CM0009618



Introduction

RPL26 Rabbit mAb is a recombinant monoclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. Recognizes human, mouse, rat RPL26 (60S ribosomal protein L26), a core large ribosomal subunit component.

[Learn More](#)

Parameter	Value
Product Name	RPL26 Rabbit mAb
Remarks/Alias	L26; uL24; DBA11
Species	Human
GeneID (Human)	6154
GeneID	6154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human RPL26 (NP_000978.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61254
Protein Weight	17kDa
Shipping	Ice bag