


Anti-Rpl26 Rabbit Polyclonal Antibody (Rpl26/Ui24) For Wb, If/Icc, Elisa - P61254

Item Number: CM0022517



Introduction

This rabbit polyclonal antibody targets human RPL26 (large ribosomal subunit protein uL24), validated for WB, IF/ICC, and ELISA. It cross-reacts with mouse and rat, making it a useful tool for ribosomal and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPL26 Rabbit pAb
Remarks/Alias	L26; uL24; DBA11; RPL26
Species	Human
Gene ID	6154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-95 of human RPL26 (NP_000978.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61254
Protein Weight	17kDa
Shipping	Ice bag