

Anti-Rpl26L1 Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q9Unx3

Item Number: CM0010495

Introduction

Rabbit monoclonal antibody against human RPL26L1 (uL24-like), detecting a 17 kDa ribosomal protein. Validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for ribosome and translation research.

[Learn More](#)

Parameter	Value
Product Name	RPL26L1 Rabbit mAb
Remarks/Alias	RPL26P1; RPL26L1
Species	Human
GeneID (Human)	51121
GeneID	51121
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RPL26L1 (Q9UNX3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNX3
Protein Weight	17kDa
Shipping	Ice bag