

Anti-Rpl26L1 Rabbit Polyclonal Antibody For Wb & Elisa - Q9Unx3

Item Number: CM0030884

Introduction

Rabbit polyclonal antibody targeting human RPL26L1 (ribosomal protein uL24-like), reactive in WB and ELISA applications. Ideal for studying ribosome structure, protein synthesis, and uL24 function in large subunit biology.

[Learn More](#)

Parameter	Value
Product Name	RPL26L1 Rabbit pAb
Remarks/Alias	RPL26P1; RPL26L1
Species	Human
Gene ID (Human)	51121
Gene ID	51121
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1-145 of human RPL26L1 (NP_057177.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UNX3
Protein Weight	17kDa
Shipping	Ice bag