


Anti-Rpl27 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P61353

Item Number: CM0013039



Introduction

Rabbit polyclonal antibody against human RPL27 (60S ribosomal protein L27, eL27), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for ribosomal processing and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPL27 Rabbit pAb
Remarks/Alias	L27; eL27; DBA16; RPL27
Species	Human
Gene ID (Human)	6155
Gene ID	6155
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-136 of human RPL27 (NP_000979.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61353
Protein Weight	16kDa
Shipping	Ice bag