


Anti-Rpl3 Rabbit Polyclonal Antibody For Wb, Elisa - P39023

Item Number: CM0026292



Introduction

Anti-RPL3 Rabbit Polyclonal Antibody for WB and ELISA. Recognizes human and mouse large ribosomal subunit protein uL3 (RPL3, TARBP-B). Ideal for ribosome assembly and translation studies. P39023.

[Learn More](#)

Parameter	Value
Product Name	RPL3 Rabbit pAb
Remarks/Alias	L3; uL3; ASC-1; TARBP-B; RPL3
Species	Human
Gene ID (Human)	6122
Gene ID	6122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RPL3 (NP_000958.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P39023
Protein Weight	46kDa
Shipping	Ice bag