

## Anti-Rpl30 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62888

Item Number: CM0017188



### Introduction

High-quality rabbit polyclonal antibody targeting human ribosomal protein L30 (RPL30). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for ribosome and translation studies.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | RPL30 Rabbit pAb                                                                                                 |
| Remarks/Alias    | L30; eL30; RPL30                                                                                                 |
| Species          | Human                                                                                                            |
| GeneID (human)   | 6156                                                                                                             |
| GeneID           | 6156                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human RPL30 (NP_000980.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P62888                                                                                                           |
| Protein Weight   | 13 kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                          |