

Anti-Rpl32 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62910

Item Number: CM0013836



Introduction

RPL32 rabbit polyclonal antibody validated for Western blot, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat RPL32. Ideal for ribosomal protein and translation studies.

[Learn More](#)

Parameter	Value
Product Name	RPL32 Rabbit pAb
Remarks/Alias	L32; eL32; PP9932; RPL32
Species	Human
GeneID (Human)	6161
GeneID	6161
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human RPL32 (NP_000985.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62910
Protein Weight	16kDa
Shipping	Ice bag