

Anti-Rpl34 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P49207

Item Number: CM0021928



Introduction

Rabbit polyclonal antibody targeting human ribosomal protein L34 (RPL34), designed for WB, IF/ICC, and ELISA applications. Detects human and mouse eL34 (13 kDa). Ideal for ribosome and translation research.

[Learn More](#)

Parameter	Value
Product Name	RPL34 Rabbit pAb
Remarks/Alias	L34; eL34; RPL34
Species	Human
GeneID (Human)	6164
GeneID	6164
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of human RPL34 (NP_000986.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49207
Protein Weight	13kDa
Shipping	Ice bag