


Anti-Rpl35 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42766

Item Number: CM0020162



Introduction

Rabbit polyclonal antibody against human RPL35 (60S ribosomal protein L35), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for ribosome research and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPL35 Rabbit pAb
Remarks/Alias	L35; uL29; DBA19; RPL35
Species	Human
GeneID (Human)	11224
GeneID	11224
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-120 of human RPL35 (NP_009140.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42766
Protein Weight	15kDa
Shipping	Ice bag