


Anti-Rpl36 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y3U8

Item Number: CM0029544



Introduction

RPL36 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat RPL36, a 60S ribosomal protein eL36. For ribosome and translation studies.

[Learn More](#)

Parameter	Value
Product Name	RPL36 Rabbit pAb
Remarks/Alias	L36; eL36; RPL36
Species	Human
GeneID (Human)	25873
GeneID	25873
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human RPL36 (NP_056229.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3U8
Protein Weight	12kDa
Shipping	Ice bag