


Anti-Rpl36A Polyclonal Antibody For Wb, Elisa - P83881

Item Number: CM0016217



Introduction

Rabbit polyclonal antibody targeting human RPL36A (eL42, L44L). Validated for Western blot and ELISA. Suitable for ribosomal protein analysis, cell growth, and migration research.

[Learn More](#)

Parameter	Value
Product Name	RPL36A Rabbit pAb
Remarks/Alias	L36A; L44L; MIG6; eL42; RPL44; RPL36A
Species	Human
Gene ID (Human)	6173
Gene ID	6173
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 37-142 of human RPL36A (NP_066357.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P83881
Protein Weight	12kDa
Shipping	Ice bag