


Anti-Rpl37 Polyclonal Antibody For Wb, Elisa - P61927

Item Number: CM0019248



Introduction

Rabbit polyclonal antibody specifically recognizing human 60S ribosomal protein L37 (RPL37), validated for Western blotting and ELISA. Cross-reacts with human, mouse, and rat. Suitable for ribosome and protein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	RPL37 Rabbit pAb
Remarks/Alias	L37; eL37; RPL37
Species	Human
Gene ID	6167
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-97 of human RPL37 (NP_000988.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61927
Protein Weight	11kDa
Shipping	Ice bag