


## Anti-Rpl39 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62891

Item Number: CM0020427



### Introduction

High-quality rabbit polyclonal antibody targeting human RPL39, validated for WB, IF/ICC, ELISA. Specific detection of the 6kDa 60S ribosomal protein L39, ideal for ribosomal research and translation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | RPL39 Rabbit pAb                                                                                                        |
| Remarks/Alias    | L39; eL39; RPL39P42; RPL39_23_1806; RPL39                                                                               |
| Species          | Human                                                                                                                   |
| GeneID (Human)   | 6170                                                                                                                    |
| GeneID           | 6170                                                                                                                    |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-51 of human RPL39 (NP_000991.1). |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                                                       |
| Cross Reactivity | Human                                                                                                                   |
| SWISS            | P62891                                                                                                                  |
| Protein Weight   | 6kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                 |