


Anti-Rpl4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P36578

Item Number: CM0027633



Introduction

High-quality rabbit polyclonal antibody against human RPL4 (uL4), validated for WB, IHC-P, and ELISA. Suitable for ribosomal protein studies across human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	RPL4 Rabbit pAb
Remarks/Alias	L4; uL4; RPL4
Species	Human
GeneID (Human)	6124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RPL4 (NP_000959.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36578
Protein Weight	48kDa
Shipping	Ice bag