


Anti-Rpl5 Polyclonal Antibody For Wb, Ihc-P, Elisa - P46777

Item Number: CM0020676



Introduction

RPL5 Rabbit Polyclonal Antibody validated for Western Blot, Immunohistochemistry-Paraffin, and ELISA in human, mouse, and rat samples. Ideal for ribosomal protein L5/uL18 research.

[Learn More](#)

Parameter	Value
Product Name	RPL5 Rabbit pAb
Remarks/Alias	L5; uL18; MSTP030; PPP1R135; RPL5
Species	Human
GeneID	6125
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human RPL5 (NP_000960.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46777
Protein Weight	34kDa
Shipping	Ice bag