


Anti-Rpl6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q02878

Item Number: CM0015049



Introduction

RPL6 Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA, reactive with human, mouse, rat. Component of the large ribosomal subunit, involved in protein synthesis and HTLV-I Tax-responsive element binding.

[Learn More](#)

Parameter	Value
Product Name	RPL6 Rabbit pAb
Remarks/Alias	L6; eL6; TXREB1; SHUJUN-2; TAXREB107; RPL6
Species	Human
GeneID (Human)	6128
GeneID	6128
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 188-288 of human RPL6 (NP_000961.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02878
Protein Weight	33kDa
Shipping	Ice bag