


## Anti-Rpl7A Monoclonal Antibody For Wb, Elisa - P62424

Item Number: CM0002991



### Introduction

Rabbit monoclonal antibody targeting human RPL7A (eL8, SURF3), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for ribosomal protein studies in cytoplasm.

[Learn More](#)

| Parameter        | Value                                                                                                        |
|------------------|--------------------------------------------------------------------------------------------------------------|
| Product Name     | RPL7A Rabbit mAb                                                                                             |
| Remarks/Alias    | L7A; eL8; TRUP; SURF3; PLA-X polypeptide; RPL7A                                                              |
| Species          | Human                                                                                                        |
| GeneID (Human)   | 6130                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human RPL7A (P62424). |
| Source           | Rabbit                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                        |
| Application      | WB, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                            |
| SWISS            | P62424                                                                                                       |
| Protein Weight   | 30 kDa                                                                                                       |
| Shipping         | Ice bag                                                                                                      |