


Anti-Rpl7A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62424

Item Number: CM0018195



Introduction

Rabbit polyclonal antibody against human RPL7A, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for ribosomal protein research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	RPL7A Rabbit pAb
Remarks/Alias	L7A; eL8; TRUP; SURF3; RPL7A
Species	Human
Gene ID	6130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of human RPL7A (NP_000963.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62424
Protein Weight	30kDa
Shipping	Ice bag