


Anti-Rpl7A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62424

Item Number: CM0016437



Introduction

RPL7A rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets 60S ribosomal protein L7a (eL8), a core component of the large ribosomal subunit essential for cellular protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	RPL7A Rabbit pAb
Remarks/Alias	L7A; eL8; TRUP; SURF3; RPL7A
Species	Human
Gene ID (Human)	6130
Gene ID	6130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of human RPL7A (NP_000963.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62424
Protein Weight	30kDa
Shipping	Ice bag