

Anti-Rpl8 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62917

Item Number: CM0016693



Introduction

Rabbit polyclonal antibody targeting RPL8 (large ribosomal subunit protein uL2). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat. Ideal for ribosome biology and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPL8 Rabbit pAb
Remarks/Alias	L8; uL2; RPL8
Species	Human
GeneID (Human)	6132
GeneID	6132
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-257 of human RPL8 (NP_000964.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62917
Protein Weight	28kDa
Shipping	Ice bag