


Anti-Rpl9 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P32969

Item Number: CM0027551



Introduction

Rabbit polyclonal antibody targeting human ribosomal protein L9 (RPL9), validated for WB, IF/ICC, and ELISA. Shows cross-reactivity with human and mouse. Ideal for ribosome biogenesis and protein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	RPL9 Rabbit pAb
Remarks/Alias	L9; uL6; NPC-A-16; RPL9
Species	Human
Gene ID (Human)	6133
Gene ID	6133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-192 of human RPL9 (NP_000652.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P32969
Protein Weight	22kDa
Shipping	Ice bag