


Anti-Rpl9 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P32969

Item Number: CM0014480



Introduction

Rabbit polyclonal antibody against human RPL9 (60S ribosomal protein L9/uL6) for WB, IF/ICC, ELISA. Cross-reacts with mouse. Suitable for ribosome and protein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	RPL9 Rabbit pAb
Remarks/Alias	L9; uL6; NPC-A-16; RPL9
Species	Human
Gene ID (Human)	6133
Gene ID	6133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-192 of human RPL9 (NP_000652.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P32969
Protein Weight	22kDa
Shipping	Ice bag