


Anti-Rplp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05386

Item Number: CM0027420



Introduction

Rabbit polyclonal antibody targeting the large ribosomal subunit protein P1 (RPLP1). Validated for Western blot (WB), immunofluorescence (IF/ICC), and ELISA. Suitable for human, mouse, rat, and zebrafish samples.

[Learn More](#)

Parameter	Value
Product Name	RPLP1 Rabbit pAb
Remarks/Alias	P1; LP1; RPP1; RPLP1
Species	Human
GeneID	6176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human RPLP1 (NP_000994.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Zebrafish
SWISS	P05386
Protein Weight	12kDa
Shipping	Ice bag