


Anti-Rplp2 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P05387

Item Number: CM0030622



Introduction

Rabbit polyclonal antibody against human RPLP2 (60S acidic ribosomal protein P2). Validated for IHC-P, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat. Ideal for protein synthesis and ribosomal research.

[Learn More](#)

Parameter	Value
Product Name	RPLP2 Rabbit pAb
Remarks/Alias	P2, LP2, RPP2, D11S2243E, RPLP2
Species	Human
GeneID (human)	6181
GeneID	6181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human RPLP2 (NP_000995.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05387
Protein Weight	12kDa
Shipping	Ice bag