

Anti-Rpn1 Rabbit Polyclonal Antibody For Wb, Elisa - P04843

Item Number: CM0019892



Introduction

Rabbit polyclonal antibody against human RPN1 (ribophorin I), a subunit of the oligosaccharyl transferase complex, validated for Western blot and ELISA with cross-reactivity to mouse and rat. Ideal for studying N-glycosylation and ER biology.

[Learn More](#)

Parameter	Value
Product Name	RPN1 Rabbit pAb
Remarks/Alias	RPN1; OST1; RBPH1; ribophorin I
Species	Human
GeneID	6184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-350 of human RPN1 (NP_002941.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04843
Protein Weight	68kDa
Shipping	Ice bag