


Anti-Rpn2 Polyclonal Antibody For Wb, Elisa - P04844

Item Number: CM0029265



Introduction

Rabbit polyclonal antibody targeting human RPN2 (ribophorin II), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for investigating N-glycosylation and ER-associated protein processing.

[Learn More](#)

Parameter	Value
Product Name	RPN2 Rabbit pAb
Remarks/Alias	SWP1, RPNII, RIBIIR, RPN-II, RPN2
Species	Human
GeneID	6185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-540 of human RPN2 (NP_002942.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04844
Protein Weight	69kDa
Shipping	Ice bag