

Anti-Rpp25 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bul9

Item Number: CM0029737



Introduction

Rabbit polyclonal antibody targeting human RPP25, a ribonuclease P/MRP subunit involved in tRNA and rRNA processing. Validated for WB, IHC-P, ELISA with human and mouse cross-reactivity. Ideal for nucleolar protein research.

[Learn More](#)

Parameter	Value
Product Name	RPP25 Rabbit pAb
Remarks/Alias	RPP25
Species	Human
GeneID (Human)	54913
GeneID	54913
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human RPP25 (NP_060263.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BUL9
Protein Weight	21kDa
Shipping	Ice bag