

Anti-Rpp38 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P78345

Item Number: CM0023220



Introduction

High-quality rabbit polyclonal antibody targeting human RPP38, validated for WB, IHC-P, IF/ICC, and ELISA. Ideal for studying ribonuclease P/MRP complexes, tRNA and rRNA processing, and nucleolar biology.

[Learn More](#)

Parameter	Value
Product Name	RPP38 Rabbit pAb
Remarks/Alias	RPP38
Species	Human
Gene ID (Human)	10557
Gene ID	10557
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human RPP38 (NP_892117.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P78345
Protein Weight	32kDa
Shipping	Ice bag