

Anti-Rpp40 Polyclonal Antibody For Wb, Elisa - 075818

Item Number: CM0019835



Introduction

Rabbit polyclonal antibody targeting human RPP40 (RNase P subunit p40). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for tRNA processing and nucleolus studies.

[Learn More](#)

Parameter	Value
Product Name	RPP40 Rabbit pAb
Remarks/Alias	RNASEP1; bA428j1.3; RPP40
Species	Human
GeneID (Human)	10799
GeneID	10799
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-210 of human RPP40 (NP_001273061)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75818
Protein Weight	42kDa
Shipping	Ice bag