


Anti-Rprd1A Polyclonal Antibody For Wb, Elisa - Q96P16

Item Number: CM0023932



Introduction

Rabbit polyclonal antibody targeting human RPRD1A, suitable for WB and ELISA. Cross-reacts with human, mouse, rat. Useful for studying transcription regulation and cell cycle control. Molecular weight ~36 kDa.

[Learn More](#)

Parameter	Value
Product Name	RPRD1A Rabbit pAb
Remarks/Alias	P15RS; HsT3101; RPRD1A
Species	Human
Gene ID (Human)	55197
Gene ID	55197
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96P16
Protein Weight	36kDa
Shipping	Ice bag