

Anti-Rprd1B Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nqg5

Item Number: CM0018554



Introduction

Rabbit polyclonal antibody targeting RPRD1B (CREPT, Kub5-Hera), suitable for WB, IHC-P, ELISA with human, mouse, rat reactivity. Involved in transcriptional regulation and tumorigenesis, ideal for oncology research.

[Learn More](#)

Parameter	Value
Product Name	RPRD1B Rabbit pAb
Remarks/Alias	K-H; CREPT; NET60; C20orf77; Kub5-Hera; dj1057B20.2; RPRD1B
Species	Human
GeneID (Human)	58490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-326 of human rprd1b (Q9NQG5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQG5
Protein Weight	37kDa
Shipping	Ice bag