


Anti-Rps10 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P46783

Item Number: CM0025862



Introduction

Rabbit polyclonal antibody against human RPS10 (40S ribosomal protein S10) validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Suitable for studies of ribosome function, translation, and Diamond-Blackfan anemia research.

[Learn More](#)

Parameter	Value
Product Name	RPS10 Rabbit pAb
Remarks/Alias	S10; DBA9; eS10; RPS10
Species	Human
Gene ID (Human)	6204
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human RPS10 (NP_001190174.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46783
Protein Weight	19kDa
Shipping	Ice bag