


Anti-Rps12 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P25398

Item Number: CM0026669



Introduction

RPS12 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat samples. Recognizes 40S ribosomal protein S12 (eS12), a key component of the small ribosomal subunit involved in ribosome biogenesis and translation.

[Learn More](#)

Parameter	Value
Product Name	RPS12 Rabbit pAb
Remarks/Alias	S12; eS12; RPS12
Species	Human
GeneID (Human)	6206
GeneID	6206
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human RPS12 (NP_001007.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25398
Protein Weight	15kDa
Shipping	Ice bag