


Anti-Rps13 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62277

Item Number: CM0017451



Introduction

Rabbit polyclonal antibody against human RPS13 (ribosomal protein S13, uS15). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for ribosome biogenesis and translation studies.

[Learn More](#)

Parameter	Value
Product Name	RPS13 Rabbit pAb
Remarks/Alias	S13; uS15; RPS13
Species	Human
GeneID (Human)	6207
GeneID	6207
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 79-151 of human RPS13 (NP_001008.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62277
Protein Weight	17kDa
Shipping	Ice bag