


Anti-Rps14 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62263

Item Number: CM0025650



Introduction

Rabbit polyclonal antibody against human RPS14 (40S ribosomal protein S14), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for ribosome biogenesis and nucleolus studies.

[Learn More](#)

Parameter	Value
Product Name	RPS14 Rabbit pAb
Remarks/Alias	S14; EMTB; uS11; RPS14
Species	Human
GeneID (Human)	6208
GeneID	6208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-151 of human RPS14 (NP_005608.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62263
Protein Weight	16kDa
Shipping	Ice bag