


Anti-Rps17 Rabbit Polyclonal Antibody For Wb, Elisa - P08708

Item Number: CM0024343



Introduction

Anti-RPS17 rabbit polyclonal antibody for WB and ELISA, validated in human, cross-reacts with mouse and rat. Recognizes 40S ribosomal protein S17, a component of the small ribosomal subunit with roles in ribosome biogenesis and Diamond-Blackfan anemia.

[Learn More](#)

Parameter	Value
Product Name	RPS17 Rabbit pAb
Remarks/Alias	S17; DBA4; eS17; RPS17L; RPS17L1; RPS17L2; RPS17
Species	Human
GeneID (human)	6218
GeneID	6218
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RPS17 (NP_001012.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P08708
Protein Weight	16kDa
Shipping	Ice bag