

Anti-Rps17 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08708

Item Number: CM0022231



Introduction

Rabbit polyclonal antibody against RPS17, a 40S ribosomal protein. Validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for ribosome and protein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	RPS17 Rabbit pAb
Remarks/Alias	S17; DBA4; eS17; RPS17L; RPS17L1; RPS17L2; RPS17
Species	Human
GeneID	6218
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RPS17 (NP_001012.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08708
Protein Weight	16kDa
Shipping	Ice bag