


Anti-Rps19 Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Elisa - P39019

Item Number: CM0011809



Introduction

High-quality rabbit monoclonal antibody targeting human RPS19 (40S ribosomal protein S19). Validated for WB, IP, IF, ELISA; cross-reacts with human, mouse, rat. Ideal for ribosome biogenesis and Diamond-Blackfan anemia research.

[Learn More](#)

Parameter	Value
Product Name	RPS19 Rabbit mAb
Remarks/Alias	DBA; S19; DBA1; eS19; LOH19CR1; RPS19
Species	Human
GeneID	6223
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 66-145 of human RPS19 (P39019).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P39019
Protein Weight	16kDa
Shipping	Ice bag