


Anti-Rps19 Polyclonal Antibody For Wb, Elisa - P39019

Item Number: CM0021727



Introduction

RPS19 Rabbit pAb is a polyclonal antibody against human 40S ribosomal protein S19, validated for WB and ELISA. Ideal for studying ribosome biogenesis and Diamond-Blackfan anemia.

[Learn More](#)

Parameter	Value
Product Name	RPS19 Rabbit pAb
Remarks/Alias	DBA; S19; DBA1; eS19; LOH19CR1; RPS19
Species	Human
GeneID (Human)	6223
GeneID	6223
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human RPS19 (NP_001013.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P39019
Protein Weight	16kDa
Shipping	Ice bag