


## Anti-Rps19Bp1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q86Wx3

Item Number: CM0014928



### Introduction

Rabbit polyclonal antibody against RPS19BP1 (AROS) for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for ribosomal biogenesis, SIRT1 regulation, and p53 signaling studies.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | RPS19BP1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | AROS; S19BP; RPS19BP1                                                                                               |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 91582                                                                                                               |
| Gene ID          | 91582                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-136 of human RPS19BP1 (NP_919307.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q86WX3                                                                                                              |
| Protein Weight   | 15kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |