

Anti-Rps20 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P60866

Item Number: CM0012047



Introduction

Rabbit monoclonal antibody against human RPS20 (40S ribosomal protein S20). Validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Ideal for ribosomal research and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPS20 Rabbit mAb
Remarks/Alias	S20; uS10; RPS20
Species	Human
Gene ID (Human)	6224
Gene ID	6224
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RPS20 (P60866)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60866
Protein Weight	13kDa
Shipping	Ice bag