


## Anti-Rps21 Rabbit Polyclonal Antibody For Wb, Elisa - P63220

Item Number: CM0020009



### Introduction

RPS21 Rabbit pAb is a polyclonal antibody targeting human RPS21, a 40S ribosomal protein essential for translation. Validated for Western blot (WB) and ELISA, cross-reacts with mouse. Excellent for studying ribosome function.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | RPS21 Rabbit pAb                                                                                                 |
| Remarks/Alias    | S21; HLDF; eS21; RPS21                                                                                           |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 6227                                                                                                             |
| Gene ID          | 6227                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-81 of human RPS21 (NP_001015.1). |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | P63220                                                                                                           |
| Protein Weight   | 9kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                          |