

Anti-Rps23 Polyclonal Antibody For Wb, Elisa - P62266

Item Number: CM0019434



Introduction

Rabbit polyclonal antibody against human RPS23 (16kDa), validated for WB and ELISA. Ideal for diagnostic manufacturers and research labs studying ribosome function and translation. SWISS: P62266.

[Learn More](#)

Parameter	Value
Product Name	RPS23 Rabbit pAb
Remarks/Alias	S23; BTDD; uS12; MABAS; MCINS; PAMAS; RPS23
Species	Human
GeneID (Human)	6228
GeneID	6228
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-120 of human RPS23 (NP_001016.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P62266
Protein Weight	16kDa
Shipping	Ice bag