

## Anti-Rps24 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P62847

Item Number: CM0017962



### Introduction

CamelBio's Anti-RPS24 Rabbit Polyclonal Antibody is validated for WB, IHC-P, IP, ELISA. Targets human 40S ribosomal protein S24 (eS24), cross-reactive with mouse and rat. Ideal for ribosome biogenesis, translation, and DBA research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | RPS24 Rabbit pAb                                                                                                  |
| Remarks/Alias    | S24; DBA3; eS24; RPS24                                                                                            |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 6229                                                                                                              |
| GeneID           | 6229                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RPS24 (NP_148982.1). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, IP, ELISA                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | P62847                                                                                                            |
| Protein Weight   | 15kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |