


## Anti-Rps25 Polyclonal Antibody For Wb, If/Icc, Elisa - P62851

Item Number: CM0015982



### Introduction

Rabbit polyclonal antibody targeting ribosomal protein S25 (RPS25/eS25). Validated for WB, IF/ICC, and ELISA applications; cross-reacts with human, mouse, and rat. Ideal for studying ribosomal protein synthesis and small subunit biology.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | RPS25 Rabbit pAb                                                                                                     |
| Remarks/Alias    | S25; eS25; RPS25                                                                                                     |
| Species          | Human                                                                                                                |
| GeneID (Human)   | 6230                                                                                                                 |
| GeneID           | 6230                                                                                                                 |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human RPS25 (NP_001019.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IF/ICC, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | P62851                                                                                                               |
| Protein Weight   | 14kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |