


Anti-Rps26 Polyclonal Antibody For Wb, Elisa - P62854

Item Number: CM0014897



Introduction

Rabbit polyclonal antibody targeting human RPS26 (40S ribosomal protein S26, 13kDa). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for ribosome and translation research.

[Learn More](#)

Parameter	Specification
Product Name	RPS26 Rabbit pAb
Remarks/Alias	S26; eS26; DBA10; RPS26
Species	Human
Gene ID (Human)	6231
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human RPS26 (NP_001020.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62854
Protein Weight	13kDa
Shipping	Ice bag