


Anti-Rps27 Rabbit Polyclonal Antibody For Wb, Elisa - P42677

Item Number: CM0026631



Introduction

Rabbit polyclonal antibody against RPS27 (40S ribosomal protein S27). Validated for WB and ELISA; reacts with human and rat samples. Target weight ~9 kDa. Useful for research on ribosome biogenesis and Diamond-Blackfan anemia.

[Learn More](#)

Parameter	Value
Product Name	RPS27 Rabbit pAb
Remarks/Alias	S27; MPS1; eS27; DBA17; MPS-1; RPS27
Species	Human
GeneID (Human)	6232
Immunogen	Recombinant fusion protein corresponding to amino acids 1-84 of human RPS27 (NP_001021.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P42677
Protein Weight	9kDa
Shipping	Ice bag