


Anti-Rps27A Polyclonal Antibody For Wb, Ihc-P, Elisa - P62979

Item Number: CM0013727



Introduction

Rabbit polyclonal antibody against RPS27A, a ubiquitin-ribosomal protein fusion involved in protein degradation and translation. Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for studies on ubiquitin signaling, proteasome pathway, and ribosome biogenesis.

[Learn More](#)

Parameter	Value
Product Name	RPS27A Rabbit pAb
Remarks/Alias	UBC; S27A; eS31; CEP80; UBA80; HEL112; UBCEP1; UBCEP80; RPS27A
Species	Human
GeneID (human)	6233
GeneID	6233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human RPS27A (NP_001170884.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62979
Protein Weight	18kDa
Shipping	Ice bag